


**PACKING LIST**

- 1 - Poly Guard II Shield
- 1 - Cap & Screen
- 1 - 2.5" to 1.5" Hose Adapter



**ASSEMBLY INSTRUCTIONS**

1. Remove the Poly Guard II Shield Assembly and all associated components from its packaging.
2. If using the 2.5" diameter hose, simply connect one end of the hose to the port on the back of the Poly Guard II Shield and the other end to either your Vaniman Large Accumulator or directly to your Vaniman Dust Collector.
3. If using the 1.5" diameter hose, install the supplied white hose adapter into the port on the back of the Poly Guard II Shield. Next, insert one end of the 1.5" diameter hose into that white hose adapter until it is secure, and connect the other side to either your Vaniman Large Accumulator or directly to your Vaniman Dust Collector.

<b>2 1/2" White Hose (VMC-1023)</b>	<b>1 1/2" Grey Hose (1450) w/ VMC-1009 end</b>
<b>Poly Guard II VMC-A398:</b>	
<p>None</p>	<p>VMC-1006A *Comes w/ the Polyguard</p>

**RECOMMENDED USE**

The Poly Guard II Shield is a dust collection accessory designed to capture and contain debris and material while you work. The container is made from durable Plexiglass with open access on both sides that allows for easy use and a mesh screen to protect from unwanted items being vacuumed. If only using for one station, it is recommended that you use either the *Vanguard 1X* or *Van-I-Vac* dust collectors. If using 2 stations, it is recommended that you use either the *Vanguard Gold 2X*, *Vanguard Platinum 2X*, or *Voyager* dust collectors for ideal suction.

**REPLACEMENT PARTS**

- Poly Guard II Replacement Shield: **VMC-1046**
- Cap & Screen Replacement: **96067**
- 2.5" to 1.5" Hose Adapter Replacement: **VMC-1006A**