

## Waxing


### Vector

#### "Sensor Touch" Programmable Waxing System

Five Programs, including the "Sensor Touch" program for precision waxing. Heats up and cools down with the touch of your finger.

Programs 1 - 4 can be pre-set to provide customized temperature control for production waxing.

- Rapid heat up
- Rapid cool down
- Various tips available
- Small PKT tip comes with each unit
- Digital microprocessor electronics
- Temperature sensing for accurate control
- Automatic selectable temperatures
- Ergonomic handpiece
- Operates in °F or °C



- Hollenback 97803
- Beavertail 97804
- Small #7 97805
- Large PKT 97802
- Small PKT 97801
- Needle 97800
- Denture Spoon 97812

Vector Electronic Waxer, Includes 1-Small PKT Tip 10150 ⚡

Extra Spatula Assembly (for "hot" tip quick change) 97809

Extra Spatula Assembly with tip combination (97809 + tip catalog number)

## Accessories

High quality stainless steel casting rings and color-coded urethane sprue formers.

### Easy Rings



| Size                 | Part Number | Size                | Part Number |
|----------------------|-------------|---------------------|-------------|
| 1 1/4" D x 1 5/8" H  | VMC-F125    | 2" D x 2" H         | VMC-F20     |
| 1 1/2" D x 1 3/4" H  | VMC-F15     | 2 1/2" D x 2 1/4" H | VMC-F25     |
| 1 3/4" D x 1 9/16" H | VMC-F175    | 3" D x 2 1/4" H     | VMC-F30     |

Custom heights available. Contact Vaniman or your dealer for information. Cutting charge may apply.

### Spectrum Sprue Formers



| Size   | Color       | Part Number | Size   | Color  | Part Number |
|--------|-------------|-------------|--------|--------|-------------|
| 1 1/4" | Purple      | VMC-SB125   | 2 1/2" | Yellow | VMC-SB250   |
| 1 1/2" | Dark Blue   | VMC-SB150   | 3"     | Orange | VMC-SB300   |
| 1 3/4" | Medium Blue | VMC-SB175   | 4"     | Red    | VMC-SB400   |
| 2"     | Light Blue  | VMC-SB200   |        |        |             |