

Micro-Abrasive Blasters

SandStorm Expert

*The micro-abrasive blasting system for the laboratory professional.
A high quality and economical blasting system for moderate to heavy production.
Low maintenance media flow system virtually eliminates wasteful down-time
and lowers maintenance costs.*



- Extra large, "Easy Fill" tanks hold up to 4 pounds of media
- Exclusive non-clogging flow design with automatic air/media purging system
- "Contact-Lens" magnified viewing window reduces eye strain and improves overall viewing area
- Blast proof resin cabinet is maintenance free
- Two hand pieces
- Includes 1 precision Tungsten Carbide Tip and 1 Crystal Embedded Tip
- Foot pedal control
- Adjustable air pressure
- 3 year warranty

SandStorm Expert (50-100 micron) 80101 ⚡
 SandStorm Expert (25 & 50-100 micron) 80102 ⚡
 SandStorm Expert (50-100 & 110-250 micron) 80103 ⚡

⚡ Available in 100, 220, 230, 240 Volts @ 50/60Hz