## MULTI STATION DUST COLLECTOR

SV

## A CENTRALIZED DUST COLLECTOR IDEAL FOR HIGH PRODUCTION LABS, SCHOOLS, OR LARGE FACILITIES.

## **FEATURES**

V-SIX

Clean Air HEPA Filtration Powerful Suction for up to 6 work areas at once Push Button Filter Cleaning Complimentary Suction Design Service Vaniman Smart System Ready Limitless Suction Plumbing Options Low Maintenance & Repair Rate

## **SPECIFICATIONS**

Reliable 330 CFM Motor Quiet 59 dB Sound (high) Durable. Weighs 136 lbs. 26"L × 25"W x 41"H Requires Dedicated 20A Circuit Manufacturer Warranty Made in California, USA



