## Vanguard Gold Mobile

## **Extraoral Chairside Suction**

Part Number: 10337







HEPA filtration of viruses, chairside biological aerosols, and moisture droplets.



Super Quiet Operation



Long lasting powerful motor.



Long 6ft flexible hose can be adjusted and stays in place.

Easy to move and operate.



Optional on & off foot pedal.



## **Product Specs**





L 14"xW 10.5"xH 14.5"







39 LBS

- Internal cartridge HEPA filter
  - o Filters 99.97% of particles down to 0.3 micron
- Internal particle filter bag
  - o Filters down to 1 Micron
- Rolling caster wheels allows easy movement
- Flex hose links can be added/removed & disinfected
- Optional exhaust mounted Odor filter
  - Can absorb mercury vapors from amalgam filling removal
- Replace bag filter as needed as it will be replaced most frequently
  - o check I time per 3 months.
- Replace HEPA filter as needed based on appearance
  - check 1 time per 6 months

