

BACTERIA
VIRUSES
AEROSOLS

WE CAPTURE WHAT YOU CAN'T SEE

As a trusted US manufacturer, Vaniman has protected dental professionals from airborne contaminants and hazards for over 30 years. We know aerosols. Always have. Always will. Let us be your partner for aerosol control and air management.

Made

in USA

Contact your equipment dealer
or visit us at www.vaniman.com.



VANGUARD MOBILE
Extraoral Aerosol Suction



PURE BREEZE
Room Air Purifier

VANIMAN
PRODUCTS THAT SOLVE PROBLEMS



VANGUARD MOBILE

EXTRAORAL SUCTION SYSTEM

**A CHAIRSIDE AEROSOL CAPTURE SYSTEM
WITH HEPA FILTRATION & QUIET OPERATION.**

FEATURES

- 360 Adjustable Hose Capture System
- Variable Speed Options
- Remote Foot On/Off Pedal
- 2 Stage HEPA Filtration
- EZ Access Maintenance Panel
- Rolling Casters
- Compact and Mobile
- Easy Disinfection & Maintenance

SPECIFICATIONS

- Strong 110 CFM/23 Kpa Motor
- Quiet 53 dB Sound (high)
- Durable. Weighs 36 lbs.
- Base Dimensions 14"L x 10.5"W x 14.5"H
- Standard 120V Wall Socket
- Manufacturer Warranty
- Made in California, USA



PURE BREEZE

HEPA AIR PURIFIER

**A HIGH VOLUME HEPA AIR FILTRATION SYSTEM
THE PURIFIERS DENTAL OFFICE & LAB AIR.**

FEATURES

- 4 Stage HEPA Filtration
- Adheres to CDC Dental Guidance
- Variable Air Volume Speed Adjustment
- Rolling Casters for Easy Movement
- Long Lasting & Reliable
- Creates Air Circulation
- Captures Viruses, Dust, Allergens, Spores
- 37..5 ACH in Standard Operatory Rooms
- Simple & Easy Maintenance

SPECIFICATIONS

- Powerful 550 CFM Motor
- Quiet 59 dB Sound (high)
- Durable. Weighs 81 lbs.
- 23"L x 21"W x 35"H
- Uses Standard Wall Outlet
- Manufacturer Warranty
- Made in California, USA