

VANGUARD

DENTAL DUST COLLECTOR SERIES

VANGUARD 1X

VANGUARD GOLD 2X

VANGUARD PLATINUM 2X

PRODUCT SPECS

SIZE

L 14" x W 10.5" x H 14.5"

(356 mm x 267 mm x 369 mm)




VANGUARD MODEL	REMOTE	REPLACEABLE FILTER BAG	TWO STATIONS	MULTI SPEED	AUTOPILOT HANDPIECE SYNC
1X	✓	✓			
GOLD 2X	✓	✓	✓	✓	
PLATINUM 2X	✓	✓	✓	✓	✓

Optional pre-filter cyclone Accumulator separates 95-99% of particulates @ 0.3 µm.



INTENDED MATERIAL USE	
ACRYLIC	✓
GYPSUM	✓
DIE-STONE	✓
ALLOY	✓
QUARTZ	✓
CERAMICS	✓
ZIRCONIA	✗
BLASTING MEDIA	✗



115V AC, 15A/5A
(220V/240V opt.)


55 dBA
max operation


35 lbs.
(15.9 kg)


Vaniman
warranty

POWER & USAGE



FILTER BAGS

Replace every 3 months,
or when 1/2 full

- Filters 97% of particles @ 1 µm
- Made by hand in California
- Easy to replace
- Also suitable for use in Non-Vaniman dust collectors

EASY ORDERING

SINGLE STATION	PART NUMBER
VANGUARD 1X	10310
VANGUARD 1X w/ACCUMULATOR	10311
VANGUARD GOLD 2X	10320
VANGUARD GOLD 2X w/ACCUM.	10321
VANGUARD PLATINUM 2X	10330
VANGUARD PLATINUM 2X w/ACCUM.	10331

TWO STATION	PART NUMBER
VANGUARD GOLD 2X	10345
VANGUARD PLATINUM 2X	10374
VANGUARD GOLD 2X w/ACCUM.	10322
VANGUARD PLATINUM 2X w/ACCUM.	10332

ASSEMBLED IN USA

VANIMAN.COM | 760.723.1498

VANIMAN