

SANDSTORM PROFESSIONAL

DESIGNED SPECIFICALLY FOR DENTAL PROFESSIONALS


This self-contained sandblaster/dust collector unit is an excellent compact option for orthodontists, periodontists, dentists, and other dental professionals who want to keep the workspace clean.



MODEL FEATURES

- Extra large 1000 mL (4 lbs.) media tank with various media size options available
- Dust fan and filter built in — no need for external dust collector
- Pen nozzle handpiece with precision blasting tip
- Foot pedal blasting control
- Bright LED for enhanced visibility
- Easy cleaning
- Easy maintenance
- Vaniman Warranty

The SandStorm Professional is designed for light dental clinic or laboratory use.



MODEL SPECS

Tank size	1000 mL
Media options	Sodium Bicarbonate, Aluminum Oxide, Glass Bead
Measurements	14" wide x 14.5" deep x 12" high
Weight	20 lbs.

PN: 80401

The SandStorm Professional requires compressed air (90 PSI at 1.5 CFM, recommended 2.5 CFM). External dust collector not necessary. Unit warranty is three (3) years. Inquire with Vaniman Manufacturing Co. for details.

Assembled in
California, USA



VANIMAN

vaniman.com
+1.760.723.1498
sales@vaniman.com