


# SANDSTORM 2

A durable dental sandblaster designed to handle medium to high levels of production with different media options for different applications. It's an excellent option for dental labs with numerous workstations or in-house clinical labs of all sizes.

## MODEL FEATURES

- Two media tanks with various media size options:
  - 25-50 micron
  - 50-100 micron
  - 110-250 micron
- Air regulator assembly for adjustable air pressure (10-90 PSI operating range)
- Two pen nozzles with precision blasting tips
- Foot pedal blasting control
- Bright LED for enhanced visibility
- Easy maintenance
- Vaniman Warranty

The SandStorm 2 is designed for medium to heavy dental laboratory use.



## MODEL SPECS

Tank size	1000 mL
Media options	Sodium Bicarbonate, Aluminum Oxide, Glass Bead
Measurements	16.5" wide x 16.5" deep x 10.5" high
Weight	21 lbs.

PN: 80301

The SandStorm 2 requires compressed air and an external dust collector (Vaniman products sold separately). Inquire with Vaniman Manufacturing Co. for details.

Manufactured in  
California, USA



# VANIMAN

vaniman.com  
+1.760.723.1498  
sales@vaniman.com