

VANGUARD MOBILE 2.0



EXTRAORAL SUCTION DEVICE FOR AEROSOLS

PART NUMBER 10337

BASE UNIT SPECS

SIZE 

L 14" x W 10" x H 21.5"
(356 mm x 254 mm x 546 mm)

Not including casters or 4' flexible capture hose

The Vanguard Mobile 2.0 provides an optimal design for lowering exposure to infected aerosols generated during normal clinical procedures.

- Whisper quiet
- Durable casters for portability
- Adjustable speed
- Compact, space-saving design
- Optional exhaust mounted odor filter (can absorb mercury vapors from amalgam filling removal)



115V AC, 15A/5A
(220V/240V opt.)



39 lbs.
(17.7 kg)



50-54 dBA
operation



Vaniman
warranty

Foot pedal start/stop included

HEPA FILTER

Filters 99.97% of particles down to 0.3 µm; replace every 6-12 months.

WOVEN PRE-FILTER

Filters 97% of particles down to 1 µm; replace every 3 months.

This chairside dental aerosol vacuum safely captures aerosols during AGPs with simple hands-free operation. Position the capture mouth within 6 inches of a patient's mouth for ideal aerosol capture.

Clinical studies have shown the Vanguard Mobile 2.0 to reduce aerosolized bacteria by as much as 84% and spatter bacteria load by as much as 95% when used in conjunction with HVE.

See clinical study and additional materials for more information.

Visit vaniman.com for details.

ASSEMBLED IN USA

VANIMAN.COM | 760.723.1498

VANIMAN