

PACKING LIST

- 1 - SandStorm XL Cabinet
- 1 - Instructions
- 1 - DC Power Supply

INSTALLATION & OPERATION

1. Remove the unit and all packaging material from the box.
2. Plug the DC Power Supply into an electrical outlet, then into the right side of the light assembly.

REPLACEMENT PARTS

<i>Fig. #</i>	<i>Description</i>	<i>Reorder Part Number</i>
1	Vent	2613
2	Blue Cap	1751
3	Floor Screen Grommet	2049
4	Floor Screen	2604
5	LED Light Assembly	97969
6	Screw, #8 x 1/4" FHSMS	2609
7	Window Assembly (with window, latch, and hinge)	2620
8	Window Only	2606A
9	Latch and Spacer	2025 & 2026
10	Static Cling Protector Shield	97951
11	Clamp	1413
12	XL Cuff (black cloth)	2033XL
13	Optional Gloves	2605
14	Window Gasket Kit	2030B
15	Red Cap	2608